


Explanation of Event Flag Labels in Data Recordings respirationAdj1 (NRO25)



At the start of a recording (in main code):
255 = if it appears in the recording, then you
know that EEG is connected to experimental
PC (gullveig).


Experimental Task / Thr1F

10 = ITI onset (Black Fixation Cross)
20 = Red Fixation Cross
SA = Stimulus Onset

First 200 | 210:
Detection Response Button Press

Second 200 | 210:
Confidence Response Button Press

(200 = Button 1 / Left)
(210 = Button 2 / Right)



Resting State

150 = Start + End of Black Fixation Cross
(Resting State Recording)

180 = Start + End of Break